

Product Name :
Heating Mantles

Product Code :
PHARMACEUTICAL SUK-
AS03LAB121A



Description :

Heating Mantles

Technical Specification :

Provided with joint less hand knitted heating nets of glass yarn capable of withstanding high temperature up to 400°C. Heating Mantle is made of one piece aluminium sheet and is logged with special glass wool. The panel is provided with on-off switch, energy control knob and pilot lamp. Supplied complete with cord and plug, to work on 220V AC 50Hz single phase.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265