

Product Name :
Ampoule Filling and Sealing Machine (Electric Operated)

Product Code :
PHARMACEUTICAL SUK-
AS03LAB123A



Description :

Ampoule Filling and Sealing Machine (Electric Operated)

Technical Specification :

The unit is built on S.S. Square pipe frame & totally encompassed with S.S. Coverings, in pleasing mat finish with all aesthetic consideration.

A large opening provided at bottom side for the ease of absolute cleaning

The pendant type electrical control panel is lifted from the machine for ease of operation & safety

No. Of Syringes

50 to 100

Range/Minute

1 to 10

Range

220V

Power Required



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265