

**Product Name :**  
Ampoule Filling and Sealing Machine (Electric Operated)

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB123A



**Description :**

Ampoule Filling and Sealing Machine (Electric Operated)

**Technical Specification :**

The unit is built on S.S. Square pipe frame & totally encompassed with S.S. Coverings, in pleasing mat finish with all aesthetic consideration.

A large opening provided at bottom side for the ease of absolute cleaning

The pendant type electrical control panel is lifted from the machine for ease of operation & safety

**No. Of Syringes**

**50 to 100**

**Range/Minute**

**1 to 10**

**Range**

**Power Required**



## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265