

**Product Name :**  
Bottle Filling Machine

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB127A



**Description :**

Bottle Filling Machine

**Technical Specification :**

Bottle Filling Machine  
Semi Automatic electric operated having filling cap. 50 to 300ml. Complete stainless steel body.  
- With single syringe  
- With double syringe



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265