


**Product Name :**  
Digital Bulk Density Apparatus

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB130A



**Description :**

Digital Bulk Density Apparatus

**Technical Specification :**

**Features:**

- 4-digit microprocessor based Pre-set counter for the setting of exact no of strokes.
- A common holder for both the cylinder drops down.
- Easy to lock cushioned holders are provided to hold the cylinder.
- In power failure, the non-volatile memory of the microprocessor remembers the last mode of operation and complete the test after the power supply is resumed.

**Technical Specification:**

- Cylinder Capacity: 50 MI. x 2 Nos.
- of Drops : Adjustable between 1 to 9999 drops.
- Display : Digital 12 mm Red LED
- Drop Height : 25 mm approx.
- Dimension : 350 W x 270 D x 280 mm H
- Weight : 1 kg. approx.



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265