

Product Name :
Digital Potentiometer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB131A



Description :

Digital Potentiometer

Technical Specification :

Features:

- Highly Stable & Accurate
- 4 Electrodes For Different Titrations
- 3 ½ Digit LED Display
- Rugged & Economical

Specifications:

Range 199.9mV

0 to \pm 1999 mV

Resolution

Repeatability

Accuracy digit

Input Impedance

Operating Temperature

Display 0.5" seven segment LED display with auto polarity indication

Power 10% AC, 50 Hz.

Dimensions 75mm

Weight(Approx.)



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265