

**Product Name :**  
Electronic Stimulator

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB133A



**Description :**

Electronic Stimulator

**Technical Specification :**

Available in four models:

1. Student Stimulators Transistorised : Out put 0-25  
DC in two steps, duration .5 to 1.5ms frequency 5 to 100 c/s in 11  
steps, modes, single, repetitive or an external trigger, for 220V 50Hz.
2. Student Electronic Stimulator Latest Design :
3. Stimulator General Purpose : Provides rectangular pulse of controlled amplitude 0.1 - 100V. Frequency  
0.1 - 1000 PPS & width 0.01 to 100ms.
4. Research Stimulator With R.F. Isolation unit :  
Provides rectangular pulse of controlled polarity, Amplitude 1.5mV to  
150V, frequency 0.1 to 10 KCS delay 10 micro sec to 1 sec & width 10  
micro sec to 1 sec.



## Astra Scientific

[www.astrascientific.com](http://www.astrascientific.com), **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265