

Product Name : I.r. Moisture Balance	Product Code : PHARMACEUTICAL SUK- AS03LAB135A
--	---



Description :
I.r. Moisture Balance
Technical Specification :
<p>A basic method for determination of moisture content, measure percent by weighing and drying simultaneously. Infra red radiation are reflected by a bulb and its temperature is controlled by a regulator.</p> <p>The balance scale is divided in terms of moisture content on wet weight basic.</p> <p>Range 100 % Accuracy $\pm 0.2\%$</p> <p>Dimension 20x15x10 cm</p> <p>Weight 20 Kg</p> <p>Sample Weight 5gm, 10gm, 25gm as required</p> <p>Power 220V, 50Hz</p>



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265