

Product Name :
Mechanical Stirrer

Product Code :
PHARMACEUTICAL SUK-
AS03LAB137A



Description :

Mechanical Stirrer

Technical Specification :

Features:

- Powerful stirring and mixing
- Maximum up to 120 liters stirring capacity
- Stepless speed control
- Ability to mix high viscosity liquids
- Wide range of stirring shafts
- High durability to chemicals

0 to 1500 rpm

Rotation

~~200W~~ / Carbon brush motor

~~5000 RPM~~ / 5000 RPM

~~160~~ Working Head

~~Speed control~~ Speed control

~~32kg~~ Weight

~~220V~~ 220/230 Volts, 50 Hz

220 / 230 Volts, 50 Hz



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265