

Product Name :
Ointment Cream Filling Machine

Product Code :
PHARMACEUTICAL SUK-
AS03LAB139A



Description :

Ointment Cream Filling Machine

Technical Specification :

Features:

Durable

Weather proof

Temperature resistance

Technical Specifications:

Power Supply 220V/50/60 Hz

Speed 5-20 RPM/min/head

Filling Accuracy ±1%

Pressure 0.4-0.5 MPa

Filling Range 25 ml 100-250 ml 200-500 ml



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265