

**Product Name :**  
Rotary Vacuum Film Evaporator

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB143A



**Description :**

Rotary Vacuum Film Evaporator

**Technical Specification :**

This model of evaporator is supplied complete with A.C. non-sparking geared drive unit, support stand fitted with jack arrangement of stainless steel rod for positioning glass assembly, variable speed regulator, thermostatically controlled stainless steel water bath fitted on a jack, double coiled glass condenser with three opening for vacuum and cooling, pear shape evaporating flask capacity 800ml. receiver flask capacity one litre with ball joint, feeding tube and vapor path tube, a complete working unit. All parts are power coated with epoxy / polyster hybrid, textured, Siemens grey finish.

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com