

**Product Name :**  
Sieve Shaker Hand Operated

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB144A



**Description :**

Sieve Shaker Hand Operated

**Technical Specification :**

Sieve Shaker (hand Operated)  
Can accommodate up to 7 sieves of 8" dia during a train or gears  
operated by a hand wheel. A heavy fly wheel ensured smooth operation.  
Supplied without test sieves and lid receiver.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265