

**Product Name :**  
Tray Dryer

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB146A



**Description :**

Tray Dryer

**Technical Specification :**

Different sizes Tray Dryer are available as follows but without trays. Chamber size in mm (WxDxH) approx tray capacity :-

1. 430 x 840 x 650mm 8 Trays 16"x32"x1 1/4" size
2. 430 x 840 x 910mm 12 Trays
3. 910 x 840 x 910mm 24 Trays
4. 960 x 1270 x 1450mm 48 Trays
5. 960 x 1700 x 1700mm 96 Trays



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265