

Product Name :
U.V. Inspection Cabinet Longwave

Product Code :
PHARMACEUTICAL SUK-
AS03LAB147A



Description :

U.V. Inspection Cabinet Longwave

Technical Specification :

U.V. Inspection Cabinet Longwave

1. U.V. Chromatography : Inspection cabinet for viewing and detecting separated substances on T.L.C. layers fitted with long wave (365nm) U.V. lamp for use on 220V A.C.
2. U.V. Chromatography : Inspection cabinet for viewing and detecting separating substances on T.L.C. layers fitted with two different wave length U.V. tubes for use on 220V A.C.
3. U.V. Inspection Cabinet (Short and Long Wave) Dual Wave : Having imported filters fitted with separate low pressure mercury vapour tubes for short U.V. (254nm) and long U.V. (365nm) plus separate lamp for visible light, the low pressure mercury vapour tubes in this model facilitate immediate use for inspection as you switch on.
4. Separate U.V. tube.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265