

<b>Product Name :</b> Vortex Mixer	<b>Product Code :</b> PHARMACEUTICAL SUK- AS03LAB148A
---------------------------------------	---



<b>Description :</b>
Vortex Mixer
<b>Technical Specification :</b>
<b>Features:</b>  Superior vibration resistant vortexing Spin compatible with microcentrifuge and cold rooms Built-in microprocessor for precise timing Built-in safety features for preventing slipping Flexibility with a wide range of accessories
<b>Specifications:</b>  BTM-502  Speed 2800 rpm in continuous mode 2800 rpm in touch mode  4.2mm

**Actuation type** continuous

**Included** cup & plastic adapter

**Max. Load**

**Power Supply** 220V/230V, 50 Hz

**Optional Accessories**

- Plate adapter with rubber disc
- Plain adapter foam
- Foam with 96X0.2 well
- Foam with 60X0.5 well
- Foam with 12mm well
- Foam with different well



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265