

**Product Name :**  
Vortex Mixer

**Product Code :**  
PHARMACEUTICAL SUK-  
AS03LAB148A



**Description :**

Vortex Mixer

**Technical Specification :**

**Features:**

Simple operation as the stable vortexing  
Spin compatible with microcentrifuge and cold rooms  
Built-in microcentrifuge for easy front-loading  
Flexibility with a broad range of accessories

**Specifications:**

ASTM-502

3000 rpm in continuous mode  
2800 rpm in touch mode

Ø 21mm

**Touch type** continuous

**Included** cup & plastic adapter

**Max.** Load

**Power Supply** 220V/230V/50 Hz

**Rubber Cup** Accessories

- Plate adapter with rubber disc
- Plain adapter foam
- Foam with 96X0.2 well
- Foam with 60X0.5 well
- Foam with 12mm well
- Foam with different well



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265