

Product Name :
Benedict Roth Recording Spirometer

Product Code :
PHARMACOLOGY SUK-
AS03LAB009



Description :

Benedict Roth Recording Spirometer

Technical Specification :

Benedict Roth Recording Spirometer (BMR) consists of following features:-

Six litre capacity, Spirometer has been provided with a unit of 4-speed electric recording Kymograph with gravity writing ink-pen.

Easily accessible valves, soda-lime container with screw connections in central chamber drain cocks to all tubes and containers.

Sampling cock provided for connecting the patient to Spirometer or atmosphere.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com