

Product Name :
Digital Respinometer

Product Code :
PHARMACOLOGY SUK-
AS03LAB010



Description :

Digital Respinometer

Technical Specification :

Polished Stainless Steel tank, with capacity of 7 manometers a side.

Digital temperature display with resolution of 0.1 deg. C.

Temperature control to within +/-0.5 deg. C.

Stainless Steel immersion heater.

Stirrer for uniform temperature.

Manometers mounted on graduated scale for reading pressure changes.

Flask with center well & side flask attached to manometer.

Oscillation of manometer approximately 90 rpm.

Powder coated exterior, with wheels for easy mobility.

Capacity 7 manometers, temperature range of 5 deg.

C above ambient to 60 deg. C.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265