

Product Name :
Isolated Organ Baths Digital

Product Code :
PHARMACOLOGY SUK-
AS03LAB012



Description :

Isolated Organ Baths Digital

Technical Specification :

Versatile and complete organ baths.

Bath chambers made out of thick Transparent acrylic sheets, chemically welded to make them leak-proof.

Temperature is achieved by 250 Watts heating elements and is microprocessor based digital temperature controlled.

A compact control panel has been provided.

Supplied complete with two sets of glass inner vessels, oxygen tubes, Warming coils, fixit levers and holders.

This bath is termed as 'Two-unit isolated bath'.

The unit is workable on 220V AC 50Hz single phase.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265