

Product Name :
Variable Interrupter

Product Code :
PHARMACOLOGY SUK-
AS03LAB019



Description :

Variable Interrupter

Technical Specification :

We also provide Tetanus set with Complete unit on two stands.

Mercury cup and magnet.

Other available recorders:-

Piston recorder

Out-flow (Float) recorder

Condons drop recorder

Inchleys drop recorder

Mechanical Strhomour

Jaquet chronometer

Jackson Enterograph

Oncometer kidney

Oncometer Heart

Hair Aesthesiometer.

Compass Aesthesiometer



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265