

Product Name : Autoclave, Portable, Stainless Steel, Wing Nut Type (Sterilizer Pressure Type)	Product Code : ANALYTICAL SUK-AS03LAB001
---	--



Description : Autoclave, Portable, Stainless Steel, Wing Nut Type (Sterilizer Pressure Type)
Technical Specification : Consistent/Seam weld development guarantees microorganisms free condition. The unit is made of the stainless steel shell, ring of aluminum, top of stainless steel joined with aluminum.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265