

|   |  |
|---|--|
| <b>Product Name :</b><br>Autoclave, Portable, Stainless Steel, Wing Nut Type (Sterilizer Pressure Type) | <b>Product Code :</b><br>ANALYTICAL SUK-AS03LAB001 |
|---|--|



|  |
|--|
| <b>Description :</b><br>Autoclave, Portable, Stainless Steel, Wing Nut Type (Sterilizer Pressure Type)   |
| <b>Technical Specification :</b><br>Consistent/Seam weld development guarantees microorganisms free condition. The unit is made of the stainless steel shell, ring of aluminum, top of stainless steel joined with aluminum. |



## Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265