

| | |
|--|--|
| Product Name : Autoclave, Portable, Aluminium, Pressure Cooker Type (Sterilizer Pressure Type) | Product Code : ANALYTICAL SUK-AS03LAB002 |
|--|--|



Description :

Autoclave, Portable, Aluminium, Pressure Cooker Type (Sterilizer Pressure Type)

Technical Specification :

Consistent development guarantees microscopic organisms free condition. The unit is made of Aluminum. Outfitted with weight measure, security valve, twofold wellbeing valve, steam discharge valve, vacuum discharge valve, elastic gasket, inward holder, internal tripod stand and gentle steel remain for Autoclave.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265