

Product Name :
Autoclave, Portable, Stainless Steel, Pressure Cooker Type (Sterilizer Pressure Type)

Product Code :
ANALYTICAL SUK-AS03LAB003



Description :

Autoclave, Portable, Stainless Steel, Pressure Cooker Type (Sterilizer Pressure Type)

Technical Specification :

Consistent development guarantees bacterial free condition. The unit is made of Stainless steel. Furnished with weight cloth, wellbeing valve, twofold security valve, steam discharge valve, vacuum discharge valve, elastic gasket, inward compartment, internal tripod stand and gentle steel remain for Autoclave.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265