

**Product Name :**  
Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)

**Product Code :**  
ANALYTICAL SUK-AS03LAB004



**Description :**

Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)

**Technical Specification :**

Twofold divider single chamber Autoclave. The inward chamber, external divider, top and ring are made of stainless steel.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265