

Product Name :
Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)

Product Code :
ANALYTICAL SUK-AS03LAB005



Description :

Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)

Technical Specification :

Twofold divider single chamber Autoclave. The inward chamber, external divider, top and ring are made of stainless steel.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265