

Product Name : Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)	Product Code : ANALYTICAL SUK-AS03LAB005
--	--



Description :
Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)
Technical Specification :
Twofold divider single chamber Autoclave. The inward chamber, external divider, top and ring are made of stainless steel.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265