

Product Name :
Autoclave, Vertical, Double Walled, Radial Locking

Product Code :
ANALYTICAL SUK-AS03LAB006



Description :

Autoclave, Vertical, Double Walled, Radial Locking

Technical Specification :

Double wall chamber Autoclave. The inner chamber and outer wall made of stainless steel, ring made of stainless steel/ mild steel.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265