

Product Name :
Autoclave, Triple Wall, Vertical, High Pressure, Nut Locking

Product Code :
ANALYTICAL SUK-AS03LAB008



Description :

Autoclave, Triple Wall, Vertical, High Pressure, Nut Locking

Technical Specification :

Triple Wall, twofold chamber Autoclave. The cover, ring, inward chamber, center chamber and external divider are made of stainless steel. The instrument is totally airtight. Outfitted with weight dressing, security valve, twofold wellbeing valve, steam discharge valve, elastic gasket, gentle steel legs for Autoclave.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265