

**Product Name :**  
Autoclave, Triple Wall, Vertical, High Pressure, Nut Locking

**Product Code :**  
ANALYTICAL SUK-AS03LAB008



**Description :**

Autoclave, Triple Wall, Vertical, High Pressure, Nut Locking

**Technical Specification :**

Triple Wall, twofold chamber Autoclave. The cover, ring, inward chamber, center chamber and external divider are made of stainless steel. The instrument is totally airtight. Outfitted with weight dressing, security valve, twofold wellbeing valve, steam discharge valve, elastic gasket, gentle steel legs for Autoclave.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265