

Product Name :
Autoclave, Triple Wall, Vertical, High Pressure, with Steam Jacket

Product Code :
ANALYTICAL SUK-AS03LAB009



Description :

Autoclave, Triple Wall, Vertical, High Pressure, with Steam Jacket

Technical Specification :

Triple divider, twofold chamber Autoclave. The cover, ring, internal chamber, center chamber and external divider made of stainless steel. Center chamber fills in as steam stock coat. The instrument is completely airtight. Furnished with weight dressing, security valve, twofold wellbeing valve, steam discharge valve, elastic gasket, gentle steel legs for Autoclave. The Autoclave can hold the in the middle of 5psi to 20 psi. Extra highlights like weight control switch, low water cut off, advanced clock, computerized temperature, alert/bell are given in this unit.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com