

<b>Product Name :</b> Autoclave, Triple Wall, Vertical, High Pressure, Radial Locking with Steam Jacket	<b>Product Code :</b> ANALYTICAL SUK-AS03LAB010
--	--



**Description :**

Autoclave, Triple Wall, Vertical, High Pressure, Radial Locking with Steam Jacket

**Technical Specification :**

Triple divider, twofold chamber Autoclave. The internal chamber, center chamber and external divider made of stainless steel, ring made of stainless steel/mellow steel. The cover is made of stainless steel fixed by outspread arms through a guiding wheel made of chrome covered gentle steel. The instrument is totally sealed. Center chamber fills in as a steam stock coat. Outfitted with weight dressing, wellbeing valve, twofold security valve, steam discharge valve, elastic gasket, inward holder, S.S container and gentle steel legs for Autoclave. The Autoclave can hold the in the middle of 5psi to 20 psi. Extra controls like weight control switch, low water cut off, computerized clock, advanced temperature, caution/ringer are given this unit.

## Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

