

**Product Name :**  
Autoclave, Horizontal, Cylindrical, Double Wall

**Product Code :**  
ANALYTICAL SUK-AS03LAB011



**Description :**

Autoclave, Horizontal, Cylindrical, Double Wall

**Technical Specification :**

The sterilizer depends on the rule that microbiological activity of immersed steam at lifted temperature is fast and exhaustive. All models deal with rule of descending uprooting of air which is most prudent strategy for getting disinfection. Development :- Double walled inward chamber, external cover and top made of stainless steel, outspread locking made of gentle steel chrome covered, ring M.S painted. Autoclave mounted on tubular steel outline with ground leveling screwed flages. Weight Control:- The weight control gadget, computerized clock, advanced temperature, low water cut off gadget and ringer on finish of cycle is fused in this unit. It economies on control utilization and diminishes the continuous opening steam discharge esteem and keeps the arrival of steam in room. Working Temperature and Pressure: Sterilizing Temperature: 1200 degree Celcius Sterilizing weight: 1.2 to 1.5 Kg/cm2 (15



## **Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265