

Product Name :
Autoclave, Horizontal, Cylindrical, Triple Wall

Product Code :
ANALYTICAL SUK-AS03LAB012



Description :

Autoclave, Horizontal, Cylindrical, Triple Wall

Technical Specification :

The sterilizer depends on the rule that microbiological activity of immersed steam at lifted temperature is quick and exhaustive. All models chip away at guideline of descending relocation of air which is most efficient strategy for getting sanitization. Development :- Triple walled with steam coat. Inward chamber and steam coat made of S.S. with watertight argon welding. Cover made of S.S.. Autoclave mounted on tubular steel outline with ground leveling screwed flages . Weight Control:- The weight control gadget is joined in all electrically worked sterilizer. It economies on control utilization and diminishes the regular opening steam discharge esteem and counteracts there rent of steam in room . Working Temperature and Pressure Sterilizing Temperature: 1200 degree Celcius Sterilizing weight: 1.2 to 1.5 K g/cm2 (15 psi-20 psi) Power: Suitable to work on 440V , 3ph, 50 Hz , Ac Supply

Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

