

Product Name :
Oven, Hot Air, Memmert Type, Stainless Steel, Thermostatic

Product Code :
ANALYTICAL SUK-AS03LAB013



Description :

Oven, Hot Air, Memmert Type, Stainless Steel, Thermostatic

Technical Specification :

Hot air stoves are appropriate for temperature up to 250°C for chemicals, dishes in labs, look into foundations, enterprises, doctor's facilities and R and D focuses.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265