

Product Name :
Oven, Hot Air, Memmert Type, Stainless Steel, Thermostatic with Fan

Product Code :
ANALYTICAL SUK-AS03LAB015



Description :

Oven, Hot Air, Memmert Type, Stainless Steel, Thermostatic with Fan

Technical Specification :

Hot air broilers are reasonable for temperature up to 250°C for chemicals, china in labs, look into establishments, ventures, doctor's facilities and R and D focuses. They are utilized as a part of different applications like warming, drying, disinfecting and preparing.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265