

Product Name :

Oven, Hot Air, Memmert Type, Stainless Steel, Digital Controller with Fan

Product Code :

ANALYTICAL SUK-AS03LAB017

**Description :**

Oven, Hot Air, Memmert Type, Stainless Steel, Digital Controller with Fan

Technical Specification :

Hot air stoves are reasonable for temperature up to 250°C for chemicals, dishes in labs, inquire about organizations, enterprises, clinics and R and D focuses. They are utilized as a part of different applications like warming, drying, sanitizing and heating.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265