

**Product Name :**

Oven, Hot Air, Memmert Type, Stainless Steel, PID Controller with Fan

**Product Code :**

ANALYTICAL SUK-AS03LAB019

**Description :**

Oven, Hot Air, Memmert Type, Stainless Steel, PID Controller with Fan

**Technical Specification :**

Hot air broilers are reasonable for temperature up to 250°C for chemicals, crystal in labs, explore establishments, businesses, doctor's facilities and R and D focuses. They are utilized as a part of different applications like warming, drying, cleaning and preparing.

**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265