

**Product Name :**  
Oven, Hot Air, Memmert Type, Aluminium, PID Controller with Fan

**Product Code :**  
ANALYTICAL SUK-AS03LAB020



**Description :**

Oven, Hot Air, Memmert Type, Aluminium, PID Controller with Fan

**Technical Specification :**

Hot air stoves are appropriate for temperature up to 250°C for chemicals, crystal in labs, examine organizations, businesses, doctor's facilities and R and D focuses. They are utilized as a part of different applications like warming, drying, disinfecting and heating.



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265