

Product Name : Oven, Hot Air, Memmert Type, Stainless Steel, PID with Computer Attachment	Product Code : ANALYTICAL SUK-AS03LAB023
---	--



Description :
Oven, Hot Air, Memmert Type, Stainless Steel, PID with Computer Attachment
Technical Specification :
It is furnished with P.I.D. controller with RS 232/RS 485 attachment for PC attachment to control the stove through PC. Hot air stoves are appropriate for temperature up to 250°C for chemicals, crystal in labs, investigate organizations, businesses, healing centers and R and D focuses. They are utilized as a part of different applications like warming, drying, sanitizing and preparing.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265