

Product Name :
Vacuum Oven

Product Code :
ANALYTICAL SUK-AS03LAB024



Description :

Vacuum Oven

Technical Specification :

It is a very flexible, smaller broiler with a huge valuable working region. It can be worked at environmental or lessened weight. The stove empowers warm touchy materials to be dried at low temperature. Then again, more steady examples might be dried all the more quickly at higher temperature, without fan dissemination. This stove is additionally especially helpful for controlled environment warming (e.g. curing of gums under nitrogen or decreased oxygen conditions) since the internal chamber is fixed. It might likewise be appropriate for dissolvable expulsion work gave that satisfactory fumes and ventilation are accommodated the stove and room individually.

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com