

Product Name :

Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic

Product Code :

ANALYTICAL SUK-AS03LAB025

**Description :**

Incubator, Bacteriological, Memmert Type, Stainless Steel, Thermostatic

Technical Specification :

Hatcheries are appropriate for assortment of employments like drying of slides, paraffin inserting, tissue culture work brooding of counter acting agent test, microbiological assurance, hatching of hydroxy steroids and crystallization thinks about. Hatcheries are reasonable for temperature up to 80°C with a precision of $\pm 1^\circ \text{C}$.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265