

Product Name : Incubator, Bacteriological, Memmert Type, Aluminium, Digital Temperature Controller with Fan	Product Code : ANALYTICAL SUK-AS03LAB028
-----------------------------------------------------------------------------------------------------------------------	----------------------------------------------------



Description :

Incubator, Bacteriological, Memmert Type, Aluminium, Digital Temperature Controller with Fan

Technical Specification :

Incubator are reasonable for assortment of employments like drying of slides, paraffin implanting, tissue culture work brooding of immune response test, microbiological assurance, hatching of hydroxy steroids and crystallization contemplates. Hatcheries are appropriate for temperature up to 80°C with a precision of +/- 1° C.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265