

**Product Name :**

Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Computer Attachment

**Product Code :**

ANALYTICAL SUK-AS03LAB029

**Description :**

Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Computer Attachment

**Technical Specification :**

Incubators are suitable for a variety of uses like drying of slides, paraffin embedding, tissue culture work incubation of antibody test, microbiological determination, incubation of hydroxysteroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of  $\pm 1^\circ \text{C}$ .



**Astra Scientific**

www.astrascientific.com, Email: info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265

