

Product Name : Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Computer Attachment	Product Code : ANALYTICAL SUK-AS03LAB029
---	--



Description :

Incubator, Bacteriological, Memmert Type, Stainless Steel, PID Controller with Computer Attachment

Technical Specification :

Incubators are suitable for a variety of uses like drying of slides, paraffin embedding, tissue culture work incubation of antibody test, microbiological determination, incubation of hydroxysteroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of +/- 1° C.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265

