

Product Name :
Incubator with Shaker, Stainless Steel

Product Code :
ANALYTICAL SUK-AS03LAB037



Description :

Incubator with Shaker, Stainless Steel

Technical Specification :

To accommodate 16 carafe of 250 ml. Hatchery with variable speed shaker gave computerized r.p.m. marker and advanced temperature controller cum pointer. It has a stainless steel council of size 450 x 450 x 450 mm.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265