

Product Name :
Hand Centrifuge, Heavy

Product Code :
ANALYTICAL SUK-AS03LAB046



Description :

Hand Centrifuge, Heavy

Technical Specification :

For 15ml. test tubes with swing out head and polyamide cans to house test tubes. It has cinching gadget for settling on table. It offers 1500 to 2500 cycles for every moment.



Astra Scientific

www.astrascientific.com, **Email:** info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265