

Product Name :
Centrifuge, Economy, Round Type

Product Code :
ANALYTICAL SUK-AS03LAB047



Description :

Centrifuge, Economy, Round Type

Technical Specification :

For 15ml. test tubes with swing out rotor head. It has 5 stage speed controller. Supply 230 volts, single stage 50Hz AC. Speed extend from 0 to 3600 RPM.



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265