

**Product Name :**  
Clinical Centrifuge, Economy

**Product Code :**  
ANALYTICAL SUK-AS03LAB048



**Description :**

Clinical Centrifuge, Economy

**Technical Specification :**

For 15ml. test tubes with settled edge rotor head. It has 5 stage speed controller. Supply 230 volts, single stage 50Hz AC. Speed extend from 0 to 3800 RPM.



**Astra Scientific**

www.astrascientific.com, **Email:** info@astrascientific.com

**Address:** K-88, 20th Street, Annanagar, Chennai, India – 600040 **Phone:** +91-8860605265