

Product Name :
Centrifuge, Digital RPM, Economy

Product Code :
ANALYTICAL SUK-AS03LAB049



Description :

Centrifuge, Digital RPM, Economy

Technical Specification :

Table best rotator with hearty development and acrylic cover. With step less controller, Digital RPM meter and 0-a hour clock and marker. Speed 0 to 4200 RPM. It can be fitted with various sort of rotor heads according to end-client necessities. Supply 230 volts/50Hz



Astra Scientific

www.astrascientific.com, Email: info@astrascientific.com

Address: K-88, 20th Street, Annanagar, Chennai, India – 600040 Phone: +91-8860605265